

glucose monitoring you can be sure of

Gina P.
Dexcom Warrior*

dexcom

dexcom one⁺



dexcom G6



dexcom G7




**Dexcom ONE+, Dexcom G6 and Dexcom G7
sensors are shown here at actual size**

Product availability varies by country.
Consult your local Dexcom representative to know
which Dexcom products are available in your region.



Whatever your patients' needs or health goals are, there is a real-time Dexcom Continuous Glucose Monitoring (CGM) System for **diabetes management**.


dexcom one+



smart and easy

Best for people with diabetes who need essential, automatic glucose information


dexcom G6



powerful and connected

Best for people with diabetes who need connectivity* or a full range of customisable alerts

dexcom G7



Smart devices sold separately.†

DISCOVER OUR CGM PORTFOLIO

	Dexcom ONE	Dexcom ONE+	Dexcom G6	Dexcom G7
SENSOR				
Up to 10 days wear time	✓	✓	✓	✓
Small sensor	✓	✓	✓	✓
All-in-one sensor and transmitter	-	✓	-	✓
Warmup	2 hours	30 minutes	2 hours	30 minutes
12-hour grace period	-	✓	-	✓
Water performance [†]	Water resistant	Waterproof	Water resistant	Waterproof
ALERTS				
Customisable sound and volume	✓	✓	✓	✓
High and low alerts	✓	✓	✓	✓
Urgent Low Soon alert	-	-	✓	✓
Delay 1st High alert	✓	✓	-	✓
Silent mode [‡]	-	-	-	✓
Alert Schedule	-	-	✓	✓
APPS				
Dexcom CGM System app	✓	✓	✓	✓
Dexcom Follow app [¶] – circle of support (max 10)	-	✓	✓	✓
Dexcom Clarity	✓	✓	✓	✓
COMPATIBILITY				
Smartphone iOS and Android [†]	✓	✓	✓	✓
Standard smart watch mode [#]	-	✓	✓	✓
Direct to Watch mode ^{**} – coming soon!	-	-	-	✓
Dexcom receiver	✓	✓	✓	✓
Seamless integration with a selection of insulin pumps ^{††}	-	-	✓	✓
WEAR LOCATION				
Back of arm	✓	✓	✓	✓
Abdomen	✓	✓	✓	✓
Upper buttocks	2–17 years old	2–6 years old	2–17 years old	2–6 years old

dexcom

Jean-Pierre S.
Dexcom Warrior*

LEARN MORE ABOUT THE DEXCOM
CGM SYSTEMS AVAILABLE IN YOUR
REGION AT [DEXCOM.COM/UK](https://dexcom.com/uk)

Jean-Pierre is a sponsored spokesperson of Dexcom.

THE HOME OF AUTOMATIC GLUCOSE READINGS¹⁻³ FOR ALL PEOPLE WITH DIABETES



ACCURACY

Accuracy you can
trust day in day out



CHOICE & CONVENIENCE

Flexible and ready
for the real world



DEDICATED SUPPORT

Comprehensive
assistance for users and
healthcare professionals

*Not all Dexcom Systems have the same connectivity. Please check connectivity options at dexcom.com/en-gb/integrations-and-connectivity. †For a list of compatible devices, visit dexcom.com/compatibility. ‡Dexcom ONE+ and Dexcom G7 sensor are waterproof and may be submerged under 2.4 meters of water for up to 24 hours. Dexcom G6 and Dexcom ONE are water resistant up to 2.4 meters of water for 24 hours when properly installed. §For up to 6 hours. ¶Users must have a compatible smartphone to set up alert schedules. This feature is not available on the receiver. **Internet connectivity required for data sharing. Following requires the use of the Follow app. Users should always confirm readings on the Dexcom ONE+, Dexcom G6 or Dexcom G7 app or receiver before making treatment decisions. *Simultaneous use of a compatible smartphone is required to receive glucose readings on a compatible smart watch. **Feature availability varies by country. ††To learn more about insulin pump integrations and compatibility with Dexcom CGM Systems, visit dexcom.com/integrate.

1 Lind M, et al. *JAMA*. 2017;317(4):379–87. 2 Polonsky WH, et al. *Diabetes Care*. 2017;40(6):736–41. 3 Beck RW et al. *JAMA*. 2017;317(4):371–378.

Dexcom, Dexcom Clarity, Dexcom Follow, Dexcom One, Dexcom Share, and any related logos and design marks are either registered trademarks or trademarks of Dexcom, Inc. in the United States and/or other countries. MAT-1792 V2.0